

Dialysis Cannula

Supercath™ NEO

Safety Fistula Catheter for Hemodialysis

Fully Encased Passive
Safety Mechanism

A gentle hemodialysis treatment
with plastic catheter tube



- *Designed to Provide Effective Dialysis and Preserve Vascular Access*
- *Made in Japan*


CE 0197



MEDIKIT

Dialysis Cannula

Supercath™ NEO Safety Fistula Catheter for Hemodialysis



Catalog Number	Catheter		Color Code	Catheter Length (mm)	No. of Side Holes	Packing
	O.D. (mm) / Gauge	I.D. (mm) / Gauge				
SP802-14 (25)	2.1 / 14G	1.7 / 16G	Pale Green	25	4	50pcs / Shelf Box 1000pcs / Shipper Box (20 Shelf Boxes)
SP802-14 (33)	2.1 / 14G	1.7 / 16G		33	4	
SP802-15 (25)	1.9 / 15G	1.5 / 17G	Blue-grey	25	4	
SP802-15 (33)	1.9 / 15G	1.5 / 17G		33	4	
SP802-16 (25)	1.7 / 16G	1.3 / 18G	White	25	4	

Please read the instructions for use carefully before using the device.

*Sterilized by E.B. (Electron Beam)



Medikit Co., Ltd.

1-13-2, Yushima, Bunkyo-ku, Tokyo 113-0034, Japan
TEL:+81-3-3839-0492 | FAX:+81-3-3839-0644 | www.medikit.co.jp/english/

Manufacturer:

Togo Medikit Co., Ltd.
17148-6, Aza Kamekawa, Oaza Hichiya,
Hyuga City, Miyazaki Prefecture 883-0062, Japan

Distributed in Europe by:

Medikit Europe GmbH
Gerbermühlstrasse 7
60594 Frankfurt am Main Germany
TEL: +49-69-8700-660-60 | FAX: +49-69-8700-660-89 | www.medikit-europe.de



Authorized Representative in the European Community:
Bernas Medical
55 Rue Stephenson, 75018 Paris, France